

Homogeneous Phantom

Standard test block specialized for digital mammography and digital breast Tomosynthesis.

## **Specifications**

- · Extreme flat PMMA block with a polished appearance
- Size: 320 x 260 x 45 mm
- Thickness tolerance ± 0.1 mm

## **Applications**

- Standard test block
- Daily check on image artefacts
- Check on the correspondence between AEC sensors
- Getting rid of possible ghosts after ghosting measurements
- Detector stability tests
- The uncorrected defective detector elements test
- Evaluation of Homogeneity of the reconstructed digital breast Tomosynthesis image

## Compatibility

- The fourth edition of the European Guidelines for Quality Assurance in Breast Cancer Screening and Diagnosis (EPQC)
- EUREF type test protocol
- EUREF Breast Tomosynthesis Quality Control Protocol
- 61223-2-6
- 61223-2-3

askforinfo@artinis.com



Einsteinweg 17

6662PW Elst

The Netherlands





+31 481 350 980



Specialized for digital mammography and digital breast
Tomosynthesis.

Page 1 of 1

Homogeneous Phantom Leaflet - rev. 2102.01