
## artinis

## Intra Oral QC Kit

The Artinis Intra Oral QC Kit complies with Belgium criteria of acceptability for dental X-ray machines for diagnosis purposes.

## **Elements**

- Test object with line patterns that consist of:
  - Standard: 5 lp/mm
  - Optional 5, 6, 7, 8, 10 lp/mm
- Aluminium step-phantom (2.5 mm, 3.5 mm, 4.75 mm, 7 mm, 9 mm)
- 2 pcs PMMA plates (50 x 40 x 10 mm)
- 1 pcs aluminium plate (50 x 40 x 2 mm)
- 1 pcs aluminium plate (50 x 40 x 1 mm)
- 1 pcs aluminium plate (50 x 40 x 0.5 mm)
- 2 pcs aluminium plate (50 x 40 x 0.2 mm)
- 1 pcs aluminium plate (50 x 40 x 0.1 mm)
- 1 pcs aluminium plate (25 x 40 x 0.1 mm)









askforinfo@artinis.com







phantoms.artinis.com

Page 1 of 1

Intra Oral QC Kit Leaflet - rev. 2102.01