

MAM-162

Phantom especially developed for acceptance testing of Digital Mammography according to the German DIN 6868-162 standard. "Image quality assurance in diagnostic X-ray departments; Acceptance test of X-ray installations for digital mammography". The outline of the phantom is close to the real breast shape. A version for constancy testing is also available.





Elements

- Basis phantom (40 mm) with:
 - Aluminum step wedge (14 steps from 0-5, 2 mm)
 - 2 rows of metal balls for check of missed tissue at chest wall side
- Test Insert PMMA with quadratic Marking
- Test Insert SDNR (Signal Difference to Noise Ratio)
- Test Insert HK (High Contrast)
- Structure plate (6 mm) with:
 - Saving for the different inserts
 - 2 rows of metal balls for check of missed tissue at chest wall side
 - Saving for aluminum step wedge
- \bullet Attenuation plates (1 x 4 mm, 1 x 10 mm, 3 x 20 mm) breast size shaped
- Attenuation plate 320 x 260 x 20 mm
- Carrying case

askforinfo@artinis.com



Einsteinweg 17

6662PW Elst

The Netherlands







Page 1 of 1