

Z-resolution Phantom

Measure the inter-plane resolution or Z-resolution of Digital Breast Tomosynthesis systems.


## Highlights


 PMMA plate of 320 x 260 x 5 mm with 25 embedded aluminium spheres of 1 mm precisely in the middle of the plate's depth.

## Compatibility

 PMMA plate of 320 x 260 x 5 mm with 25 embedded aluminium spheres of 1 mm precisely in the middle of the plate's depth. Measure the inter-plane resolution or z-resolution of Digital Breast Tomosynthesis systems.







askforinfo@artinis.com







phantoms.artinis.com

Page 1 of 1

Z-resolution Phantom Leaflet - rev. 2102.01